



---

## Reverse Transfections of Adherent Cells in 96-well Plates

---

Last Revision: 2017-06-23

Version: 1.0

Written by: Mehdi Bouhaddou

Validated by / Date:

This protocol is designed as a rapid alternative to standard transfection that does not require plating cells the day before transfection. Instead, a suspension of cells is added directly to liposome complexes prepared in a 96-well plate. This protocol has been used successfully with the cells and conditions outlined below. This reverse-transfection has also been successful using larger well plates, such as 6-well plates (scaled appropriately).

This protocol is based on that from Life Technologies called: “Alternate rapid protocol for 96-well transfections without pre-plating cells”

[http://tools.lifetechnologies.com/Content/SFS/ProductNotes/F\\_Lipofectamine%202000b-040923-RD-MKT-TL-HL050602.pdf](http://tools.lifetechnologies.com/Content/SFS/ProductNotes/F_Lipofectamine%202000b-040923-RD-MKT-TL-HL050602.pdf)

1. Make DNA master mix
  - a. For each plasmid to be transfected, prepare a separate, sterile 1.5 mL centrifuge tube.
  - b. Add Opti-MEM (Life Technologies Cat #31985070) to each tube. Use 25  $\mu$ L per well to be transfected. Multiply by 1.1 to ensure enough volume is left during dispensing later.
    - i. For example, if 3 wells are to be transfected with one plasmid, put 75\*.1.1  $\mu$ L of Opti-MEM into the tube.
    - ii. We have also used regular DMEM in place of Opti-MEM.
  - c. Add DNA to each tube.
    - i. The original protocol calls for 320ng of DNA per well. We found that to be too much for most plasmids. Often, 50ng is sufficient, but this should be decided case-by-case.
  - d. Aliquot 25  $\mu$ L of the DNA master mix (es) to the 96 well plate as appropriate.
2. Make Lipofectamine master mix.
  - a. As in Step 1a above, add Opti-MEM to each tube.

- b. Add 0.6  $\mu\text{L}$  Lipofectamine™ 2000 (Invitrogen Cat #11668-019) per 25 $\mu\text{L}$  Opti-MEM (i.e. per well).
          - i. For example, for transfecting 12 wells in a 96 well plate, one would add 7.2  $\mu\text{L}$  lipofectamine to 300  $\mu\text{L}$  OptiMEM. As above, we typically prepare 10% extra for contingency.
          - ii. The ratio of lipofectamine to DNA can be optimized by altering the amount of lipofectamine added (above protocol lists 0.4 to 0.8  $\mu\text{L}$  per well).
        - c. Incubate at room temperature for five minutes.
      3. Prepare transfection complexes.
        - a. Add 25  $\mu\text{L}$  of the Lipofectamine™ 2000/OptiMEM mix (from step #2) to each well containing DNA (from step #1), slowly and dropwise. Mix gently by rocking the plate with hands (like a nutator).
        - b. Incubate at room temperature for 30 min to allow DNA-Lipofectamine™ 2000 complexes to form.
          - i. We have found that 20 minutes is often sufficient.
      4. Prepare Cells.
        - a. Lift cells of interest (see Subculture protocol).
        - b. Resuspend to a density of 10,000 cells / 150  $\mu\text{L}$  in growth medium **without** antibiotics.
          - i. This number is typical for us in a 96-well plate for somewhat sparse confluency. This of course can be optimized for downstream assay and other contingencies.
          - ii. Again, 10% extra cell suspension volume is advisable, particularly with FBS as it can cause bubbles.
        - c. Gently add 150  $\mu\text{L}$  of the cell suspension to each well containing the DNA-Lipofectamine™ 2000 complexes (from step #3) dropwise. Mix gently as above.
          - i. You should have a total of ~200 $\mu\text{L}$  of liquid in each well after this step.
      5. Incubate at 37°C in a humidified 5% CO<sub>2</sub> incubator until ready to assay (24-48 h post transfection). It is not necessary to remove the complexes or change the medium. Cells will adhere as usual in the presence of the complexes.
        - a. If the risk (or cost) of contamination is high, one can add PenStrep (1X penicillin/streptomycin, 10,000 I.U./mL / 10,000  $\mu\text{g}/\text{mL}$ , Corning Cat #30-002-CI) ~4-6 hours post transfection. To do this, one can simply replace the non-PenStrep containing media with PenStrep (1X) containing media.
          - i. This will result in the removal of the transfection reagents, but 4-6 hours should be enough time for transfection to occur successfully. This time can be extended if needed.
          - ii. Most cell lines attach sufficiently to resist aspiration after 4-6 hours, but some may not.